

ANNOUNCING OUR NEW AUTOCLAVABLE CORNEAL SHIELDS

Manufacturers of Products of Exceptional Quality Without Compromise



Improve your bottom line with our new autoclavable shields which can be reused up to 20 times or more. This reduces the cost per use to \$5.00 or less. If you are operating around the eyes, "Our Shields Protect Your Patient's Eyes and You."

Concept and Design: *Your Satisfaction is Guaranteed or we will give you a full refund!*

Comfortable: The high polish finish of our autoclavable shields provides the utmost comfort and protection along with per-cycle savings. (available in two opaque colors, grey or white).

Reusable: You may autoclave the shields up to 20 times or more. However, for safety if the surface becomes rough or scratched, they must be replaced. Order the shields with or without a handle 2 to a package. The non-handle style shields include an autoclavable silicone suction cup for easy insertion and removal of the shield. With either style they are ready to autoclave and use for \$100 per pair. Save by ordering five pair for \$450 or six pair for \$500.

USE INSTRUCTIONS: Autoclave the shields and *silicone suction cup following the instructions enclosed with the package. Once sterilized and used, the shields and suction cup should be thoroughly washed with hot water and soap, then dried. The shields and suction cup must be packaged in an autoclaving pouch, or they may be damaged. Other acceptable methods of sterilizing are ETO ethylene oxide gas, or Sterrad® gas plasma. **Do Not Flash Sterilize** or the shields will be destroyed.

Caution: Corneal edema may occur if shields are left in place over two hours. If it is necessary to leave the shields in place longer than two hours, the shields should be removed and the corneas hydrated with sterile saline before reinserting.

The sizes available are:

- **Pediatric:** Designed for infants and very small children: Dimension 18x14 mm.
- **Small:** Suggested for children and small adults: Dimension 21x18 mm.
- **Medium:** Suggested for women and smaller male globes: Dimension 22x21 mm.
- **Large:** Designed for men, vaulted for cornea clearance: Dimension 24x22 mm.

Packaging:

- The shields are packaged non-sterile in a pouch ready for autoclaving.
- Each package includes instructions for cleaning, preparing for reuse and the methods recommended for sterilizing.

*** The silicone suction cup may require replacement sooner than the shields**

Kolberg Ocular Supplies, Inc.

Call toll free: 866-259-5004 ~ Fax: 866-950-4704

Email: corneashields@aol.com ~ Web: www.CorneaProtectors.com